

# NEW!! Dermal Biopsy Punch


## ELLIPTIPUNCH™



The Rotary ELLIPTIPUNCH™  
Patent pending

Inventor, Sam Huot, MD

To order visit [DermalPunch.com](http://DermalPunch.com)



Retractable, Rotating Blade

### Features:

- Elliptically shaped biopsy for better closure
- Rotating Elliptical Blade for easy incision
- Retractable Blade for safety
- Ribbed handle for improved gripping
- Skin gripping teeth on the barrel
- Ultra-sharp stainless steel blade
- Ergonomic design
- Size engraved on the unit
- Individually packaged
- Disposable

### Sizes:

- 8mm x 20mm—soon
- 10mm x 25mm
- 12mm x 30mm—soon



N50 W13740 Overview Drive, Suite A ■ Menomonee Falls, WI 53051  
Ph: 262 373 1700 ■ Toll Free: 866 212 8466 ■ Fax 262 373 1800

Creating Innovative New Instruments for the Medical Industry

[www.DermalPunch.com](http://www.DermalPunch.com)

