



Contact: Anne Ward
Tel: (262) 373-1700
Fax: (262) 373-1800
Email: Anne@huotinstruments.com

FOR IMMEDIATE RELEASE

HUOT OFFERS CUTTING EDGE BIOPSY PUNCHES

Huot Instruments Offers Products Unique Within the World of Dermal Biopsy Punches.

Huot Instruments of Menomonee Falls Wisconsin provides doctors with alternatives to standard dermal biopsy punches currently on the market. Among the available products are the VisiPunch, a Punch with viewing windows in the blade, and the ElliptiPunch, a punch with a rotating elliptical blade.

The VisiPunch is a circular dermal biopsy punch with "*windows*" in the blade that enable the medical clinician to clearly visualize the lesion being removed to ensure a higher rate of clear surgical margins. The window also allows for easy removal of any specimen which may become lodged within the punch.

The ElliptiPunch, another inventive dermal biopsy punch, has an elliptically shaped blade; it is unique in that the blade rotates within a plastic sleeve. The elliptically shaped incision is easier to close and makes the resulting scar from the biopsy less noticeable.

Huot Instruments also offers the DermaPunch for those more comfortable with tradition. This classic disposable biopsy punch is refined using a proprietary process making it the the sharpest punch on the market. The DermaPunch is Huot's punch for the more cost conscious consumer.

Within the world of Dermal Biopsy punches Huot Instruments offers some truly unique products. The VisiPunch with its windowed blades, the ElliptiPunch with its rotating elliptical blade, and the classic DermaPunch with its razor-fine edge all address important needs within the industry.

When asked why he started the business, Dr. Sam Huot said, "Although I'm a family physician, I do a lot of biopsies. I thought an elliptical shaped punch would make it easier to remove larger lesions."

Huot Instruments is committed to researching and developing innovative medical devices for professional use.

###

For more information on these products, or to schedule an interview, please call Anne Ward at (262) 373-1700 or e-mail her at Anne@huotinstruments.com.

