


Dermal Biopsy Punch



Wittenberg VISIPUNCH™
Patent pending

Inventor, Greg Wittenberg, MD

To order visit DermalPunch.com



Features:

- Viewing Windows in the blade to “See the Lesion”
- Easy sample removal through the window
- Ultra-sharp stainless steel blade
- Perfect circular biopsy incision
- Ribbed handle for better gripping
- Ergonomic design
- Disposable
- Individually packaged
- Gamma sterilized
- CE Marked

Sizes:

- 2mm
- 3mm
- 4mm
- 5mm
- 6mm
- 8mm



N50 W13740 Overview Drive, Suite A ■ Menomonee Falls, WI 53051
Ph: 262 373 1700 ■ Toll Free: 866 212 8466 ■ Fax 262 373 1800

Creating Innovative New Instruments for the Medical Industry

www.DermalPunch.com

